


## MedVision Leonardo Adult Patient Simulator

**A completely wireless adult patient simulator designed for a wide variety of clinical training scenarios, with the ability to mimic bodily functions realistically, a lifelike aspect and tactile feel, natural full articulation, anatomically correct landmarks and airways, and automatic detection and logging of student actions**

- 
- From advanced life support to critical care
  - Transport and patient handling
  - CPR, intubation and ventilation
  - Realistic musculoskeletal structure and seamless lifelike skin
  - Highly detailed realistic face with hair
  - Visually and tactilely realistic airways
  - Measuring of blood pressure and oxygen saturation
  - ECG and Defibrillation using real devices
  - Debriefing with video recording
  - Auscultation and a comprehensive library of heart, lung and bowel sounds
  - Spontaneous breathing and vocal sounds
  - Various complications: tongue edema, pharyngeal swelling, laryngospasm, tongue fallback airway obstruction, cyanosis, convulsions, and more
  - Bleeding arm wound and bleeding leg wound modules
  - Can be used with real mechanical ventilators
  - 1 day of installation, setup and tutorial
  - 1 year of warranty

### **Key features**

A wide variety of clinical training scenarios

Various pre-installed themed medical training scenarios and patient conditions

Autonomous power supply (battery)

Battery life in continuous operation: 8 hours

Wired charging from a 110-220V power outlet

Interface in English

Simulator is a life-size model of an adult

Height of 6'

Weight of 160 lbs

Realistic weight of limbs and head

Highly realistic head with detailed facial features such as blemishes, wrinkles, pores and hair on scalp and eyebrows

Proprietary silicone skin has realistic appearance and is tactilely lifelike with a natural texture, suppleness and elasticity

Simulator silicone skin is easy to clean and will not stain from moulage, makeup or washable markers

Simulator skin covers entire body, including articulations

Patient simulator has a realistic imitated musculoskeletal structure and has the following anatomical landmarks:

- Finger bones, wrists, ankles and elbows
- Ulna, radius and tibia

|  |
|--|
| <ul style="list-style-type: none"> <li>• Palpable ribs</li> <li>• Palpable shoulder blades</li> <li>• Kneecaps mobility</li> <li>• Palpable pelvic bones</li> <li>• Palpable muscle on arms, legs and torso</li> </ul>   |
| Natural full articulation of neck, limbs and lumbar  |
| <p>Neck mobility (in degrees):</p> <ul style="list-style-type: none"> <li>• Left and right tilt: <math>\pm 45^\circ</math></li> <li>• Forward tilt: <math>+ 20^\circ</math></li> <li>• Backwards tilt: <math>- 45^\circ</math></li> <li>• Head rotation left and right: <math>\pm 45^\circ</math></li> </ul>   |
| Bendable lumbar region   |
| <p>Mobility in all joints (in degrees):</p> <ul style="list-style-type: none"> <li>• Shoulders X: <math>0\dots+180^\circ</math></li> <li>• Shoulders Y: <math>0\dots-100^\circ</math></li> <li>• Elbow: X: <math>0\dots+110^\circ</math></li> <li>• Wrist: X: <math>+45\dots+80^\circ</math></li> <li>• Forearm rotation: <math>0\dots+180^\circ</math></li> <li>• Knees: <math>0\dots-90^\circ</math></li> <li>• Ankle joint: <math>\pm 15^\circ</math></li> <li>• Hip joint: <math>0\dots+90^\circ</math></li> </ul> |
| Male and Female interchangeable genitalia  |
| Catheterization procedure (automatically logged into software) and urination (variable strength)   |
| Programmable blinking rate and eyelid position (bilateral or unilateral)   |
| Pupillary light response is automatic (pupillary reflex can be disabled to mimic trauma)   |
| Programmable dilation and constriction of pupils (bilateral or unilateral)   |
| Programmable tonic and clonic convulsions  |
| Head tilt tracking (action automatically logged into software)   |
| Chin lift tracking (action automatically logged into software)   |
| Jaw thrust tracking (action automatically logged into software)  |
| Oculocardiac reflex (action automatically logged into software)  |
| <b>Sounds</b>  |
| Auscultation areas: 21   |
| Anterior lung auscultation: 5 areas  |
| Posterior lung auscultation: 6 areas   |
| Heart auscultation sites: 5 points (aortic, pulmonic, Erb's point, tricuspid and mitral)   |
| Bowel auscultation: 4 areas  |
| Korotkoff sounds auscultation: 1 point   |
| Auscultation of normal and abnormal lung sounds: bronchovesicular respiration, vesicular respiration, diminished vesicular respiration, coarse crackles, fine crackles, wheezes, rhonchi, gurgling rhonchi, stridor, pleural friction rub, amphoric respiration, harsh respiration, asthma, covid-19, right sided pneumothorax, left sided pneumothorax, pneumonia, atelectasis, chronic bronchitis, pulmonary edema, emphysema.   |


|   |
|---|
| Auscultation of normal and abnormal heart sounds: Functional murmur, diastolic murmur, opening snap, holosystolic murmur, early systolic murmur, continuous murmur, Austin Flint murmur, pericardial rub, Graham Steell murmur, aortic valve regurgitation, aortic valve stenosis, aortic stenosis and regurgitation, congenital aortic stenosis, mitral valve regurgitation, mitral valve stenosis, mitral valve prolapse, mitral stenosis and regurgitation, pulmonary valve stenosis, pulmonary valve regurgitation, tricuspid valve regurgitation, coarctation of the aorta, hypertrophic cardiomyopathy, patent ductus arteriosus, atrial septal defect, acute myocardial infarction, congestive heart failure, systemic hypertension, acute pericarditis, dilated cardiomyopathy, pulmonary hypertension, tetralogy of Fallot, ventricular aneurysm, Ebstein's anomaly. |
| Auscultation of normal and abnormal bowel sounds: Normal bowel sound, hyperactive sounds, hypoactive sounds, borborygmus, capotement, peritoneal friction rub, normal bowel sound with bruits, irritable bowel syndrome, diarrhea, bruits due to renal artery stenosis, constipation, ulcerative colitis, Crohn's disease, paralytic ileus.   |
| Vocal sounds: pre-recorded words, phrases, coughs, screaming, moans, vomiting, difficulty breathing, and other sounds   |
| Bluescope stethoscope adapter for clean and clear auscultation sounds   |
| Use of a built-in microphone for wireless communication or patient speech   |
| <b>Cardiovascular system</b>  |
| Pulse palpation: 14 points  |
| Pulse points: carotid bilateral, brachial bilateral, radial bilateral, femoral bilateral, popliteal bilateral, dorsalis pedis bilateral, posterior tibial bilateral   |
| Pulse is activated when palpated  |
| Each pulse palpation is automatically logged into software (action log)   |
| Configurable pulse strength   |
| Defibrillation using real defibrillators (procedure details logged in software)   |
| Realistic physiological reaction when defibrillated (patient jolts)   |
| Defibrillator discharge is registered in the activity log (logged in software)  |
| ECG electrodes correct positioning tracking (logged in software)  |
| Efficient chest compressions lead to mimicked blood circulation, cardiac output, central and peripheral blood pressure and palpable pulse   |
| Cyanosis visible on hands and nasolabial triangle   |
| <b>Injections</b>   |
| IV access port, pre-installed catheter on arm   |
| IV administration with fluid or air (wet or dry)  |
| Automatic drug recognition for IV injections (logged in software)   |
| Automatic dose and injection speed recognition for IV injections (logged in software)   |
| Virtual medications: comprehensive library of medications with feedback-based pharmacokinetic response  |
| Medicine kit with syringes, drugs represented in RFID tags that attach to syringes, RFID writer (customizable)  |
| IO tibial site  |
| IO humeral head site  |
| IO administration with fluid or air (wet or dry) on humeral head site   |

|   |
|---|
| <b>Secretions</b>   |
| Forehead perspiration (diaphoresis)   |
| Lacrimation   |
| Ear secretions  |
| Nasal secretions  |
| Oral secretions   |
| Urination and catheter urination (action logged in software)  |
| Bleeding arm wound (attachable module)  |
| Bleeding leg wound (attachable module)  |
| Tourniquet application on arm and leg to stop bleeding (action logged in software)  |
| <b>Respiratory system, airway, thorax</b>   |
| Mechanical ventilator compatible (max tidal volume 800ml)   |
| Right and left lungs are fully independent  |
| Airway (mouth, oropharynx, larynx, trachea) and esophagus are based on a real patient data and anatomically correct, with lifelike tactile feel and aspect for true-to-life intubation experience |
| During spontaneous ventilation the patient simulator is breathing with automatically controlled respiratory volume and respiratory rate, keeping eucapnia and normal oxygen saturation            |
| Laryngospasm - visually lifelike and programmable   |
| Pharyngeal edema - visually lifelike and programmable   |
| Tongue edema - visually lifelike and programmable (0%/50%/100%)   |
| Trismus (lockjaw) – visually lifelike and programmable  |
| Simulated tongue fallback airway obstruction (breathing resumes with correct head positioning)  |
| Tension pneumothorax (left side or right side)  |
| When unilateral pneumothorax is active, chest excursions are accordingly unilateral, auscultation sounds are accordingly unilateral   |
| Needle decompression with audible air escape (left side or right side)  |
| Pleural cavity for chest tube insertion (bilateral)   |
| Chest compression   |
| Independent left/right airway obstruction   |
| Nasal intubation  |
| Airway intubation with an 8 mm tube (action logged in software)   |
| Endobronchial right and left lung intubation (action logged in software)  |
| Esophageal intubation with stomach distension (action logged in software)   |
| Laryngeal mask placement  |
| Cricothyrotomy and tracheostomy   |
| Consumable materials for performing cricothyrotomy and tracheostomy: neck skin, cricoid cartilage   |
| Different cricoid cartilage types: more prominent/less prominent  |
| Chest excursion when ventilated with a BVM, even when the simulator is off  |
| Software suite includes: <ul style="list-style-type: none"> <li>• Instructor software</li> </ul>  |

|  |
|--|
| <ul style="list-style-type: none"> <li>• Bedside Monitor software</li> <li>• Scenario Constructor software</li> <li>• ECG Editor software</li> <li>• Debriefing software</li> </ul>  |
| <p><b>Instructor software</b></p> <p>Allows working with pre-constructed scenarios: automatic, manual or themes; monitoring scenario progress and performance assessment; introducing various medical complications; provides detailed debriefing reports after every scenario</p> |
| Modular software structure with various add-ons  |
| Activity log for all performed actions   |
| Automatic scenario mode (a triggering action changes patient condition)  |
| Manual scenario mode (patient condition and vital signs are changed manually)  |
| Themes scenario mode (manually shift from one condition or state to another)   |
| Pre-installed themes   |
| Pre-installed simulated conditions   |
| Patient's current physiological parameters monitoring and control block  |
| Patient's bedside monitor control block  |
| Activity log for all actions performed during the session  |
| Trends display block, with forecasting graphs of how the scenario will develop   |
| CPR assessment monitoring, configuring and activating block  |
| Built-in debriefing system with video recording  |
| The debriefing system allows viewing a completed exercise video  |
| Detailed CPR statistic report  |
| CPR assessment parameters (logged in software)   |
| Chest compression depth (logged in software)   |
| Chest compression rate (logged in software)  |
| Chest compression pause length (logged in software)  |
| Time of compressions (logged in software)  |
| Ventilation rate (logged in software)  |
| Ventilation volume (logged in software)  |
| Excessive ventilation (logged in software)   |
| Saving CPR performance report to a separate file   |

|  |
|--|
| Current connection status diagram displaying block   |
| Simulator battery charge indicator   |
| Multi-language interface   |
| User database: track, create new users, edit and delete users  |
| Create, edit and delete user groups  |
| <b>Bedside Monitor software</b><br>Installed on a separate all-in-one PC – imitates a real patient monitor. Displays patient's vital signs, blood pressure monitoring, display of all ECG leads, drugs administration, virtual defibrillation, CPR monitoring. Also allows viewing additional information about the patient (MRI, CT scans, medical history) |
| Vital signs monitoring on both the bedside monitor and on the instructor's laptop  |
| Monitoring parameters presetting (anesthesiology, resuscitation, transportation, cardiology)   |
| Simultaneous waves display (3 to 7 graphs)   |
| Alarm system with the ability to mute or pause the alarm for 1 min   |
| Large digits for HR, RR, SpO2, BP, EtCO2   |
| Pulse oximeter for measuring SpO2  |
| Monitoring channels:   |
| – 12 ECG leads: I, II, III, aVL, aVR, aVF, V1, V2, V3, V4, V5, V6  |
| – Respiratory rate (RR)  |
| – Systolic blood pressure (SysBP)  |
| – Diastolic blood pressure (DiaBP)   |
| – Respiratory frequency (fR)   |
| – Respiratory curve (Resp)   |
| – Blood oxygen saturation (SpO2)   |
| – Body temperature   |
| – Non-invasive blood pressure (NIBP)   |
| – End-tidal CO2 (EtCO2)  |
| – Central venous pressure (CVP)  |
| – Pulmonary artery pressure (PAP)  |
| Panel of additional functions: blood pressure monitoring, TOF monitoring, display of all ECG leads, drugs administration, defibrillation, patient's data, CPR monitor display  |
| <b>Scenario Constructor software</b><br>Allows creation of various clinical scenarios, formation of action checklist, setting shifts and triggers between conditions using a database of events and actions, saving the script and using it with the simulator   |
| Allows creation of various proprietary training scenarios  |
| Unlimited number of states   |
| Scenario constructor is designed as a canvas for easy visualization and scenario construction  |
| Patient's physiological condition parameters settings block  |
| Number of displayed physiological parameters: 15 parameters  |

|  |
|--|
| Shifts and triggers settings block (with a database of possible actions and medication)              |
| Listing of actions and triggers required for state shift   |
| Saving patient's condition as a separate file  |
| Saving a scenario as a separate file   |
| <b>ECG Editor software</b>   |
| Allows constructing ECG graphs with their further integration into the database for use in scenarios |
| Designed as a canvas for easy visualization and graph construction                                   |
| ECG graphs pointwise construction for a certain lead   |
| Change of the interpolation between the listed points (linear, sine-wave, cubical, etc.)             |
| Saving, editing, deleting, renaming the ECG graphs for a certain lead                                |
| Impulse duration setting   |
| Real ECG representation as the graphic canvas background   |
| Received ECG curve viewing window  |
| <b>Debriefing software</b>   |
| Allows viewing the results of the exercises, discussion and analysis of the results obtained         |
| Exporting any report of a completed session  |
| Window displaying video with controls (start, stop, repeat)  |
| Bedside monitor parameters window  |
| Window displaying dynamic changes of parameters as graphs  |
| Window displaying CPR parameters as a graph  |
| Detailed CPR assessment window   |
| CPR assessment printing  |
| <b><u>Package contents:</u></b>  |
| Adult patient simulator: 1 unit  |
| Simulator battery charger: 1 unit  |
| Patient clothing (shorts)  |
| Drug injection syringe imitators: 5 syringes   |
| Holder for syringes and medication RFID tags: 1 unit   |
| Bluescope stethoscope adapter: 1 unit  |
| Bluescope stethoscope adapter charger: 1 piece   |
| Lubricant spray (for intubation)   |
| Blood Pressure Cuff: 1 unit  |
| Replaceable IV catheters: 5 pieces   |
| IO bone consumable replacement key: 1 unit   |
| IO tibia bone consumable: 10 pieces  |
| IO skin cover band: 2 extra pieces   |
| Key for replacing genitals: 1 piece  |
| IO humeral head consumables: 5 extra pieces  |



|   |
|---|
| Plugs for the urination port: 2 pieces                    |
| Pulse oximeter: 1 unit                                    |
| Neck skin: 2 extra pieces                                 |
| Chest tube insertion consumables : 6 extra pieces (3L/3R) |
| Cricothyrotomy consumables set: 2 sets                    |
| Prominent cricoid (preinstalled): 1 piece                 |
| Non-prominent cricoid: 1 piece                            |
| Refill bottles: 1 set                                     |
| Refill bottles holder: 1 piece                            |
| Leg wound module (capable of bleeding): 1 unit            |
| Arm wound module (capable of bleeding): 1 unit            |
| Wi-Fi router: 1 unit                                      |
| PC instructor: 1 unit                                     |
| PC patient monitor: 1 unit                                |
| Instructor software: 1 license                            |
| Scenario Constructor software: 1 license                  |
| Bedside Monitor software: 1 license                       |
| ECG Editor software: 1 license                            |
| Debriefing software: 1 license                            |

MedVision  
[sales@medvisiongroup.com](mailto:sales@medvisiongroup.com)  
[www.MedVisionSim.com](http://www.MedVisionSim.com)