

Cardinal Health™ Kangaroo™ Enteral Feeding Pumps

Advanced technology for reliable nutrition delivery

By your side and beyond

Supporting patients from hospital to home

The Cardinal Health™ Kangaroo™ Enteral Feeding Pump portfolio remains one of the trusted choices for patients, caregivers and clinicians looking for a dependable and versatile range of short- and long-term enteral feeding solutions. Built on a legacy of cutting-edge technology and safety, Kangaroo™ continues to shape the future of enteral feeding with the new Kangaroo OMNI™ Enteral Feeding Pump.

Enteral Feeding Pump	Kangaroo™ ePump	Kangaroo™ Joey	Kangaroo™ Connect	Kangaroo OMNI™
Feed only	✓	✓	✓	✓
Feed and flush	✓	✓	-	✓
Thick formula	-	-	-	✓
Serial port connectivity	✓	✓	✓	✓
Feed only and feed/flush accuracy	+/- 10%	+/- 7%	+/- 5%	+/- 5%
Dimensions	16.8 x 16.3 x 11.7 cm	13 x 10.4 x 9.1 cm	15.4 x 9.9 x 4 cm	17.3 x 7 x 13.5 cm
Weight	1,100 g (2.34 lbs.)	770 g (1.7 lbs.)	330 g (0.7 lbs.)	550 g (1.2 lbs.)
Operating noise (from 1 meter)	47 dB	42 dB	37 dB	37 dB
Attitude independent	✓	✓	✓	✓
Rubberized outer casing	-	✓	✓	✓
Water resistance	IPX1	IPX4	IPX6	IPX6
Display type	Monochrome LCD	Monochrome LCD	Color LCD	Color LCD
Feeding history	3 days	3 days	3 days	30 days
Battery life	15 hours	18 hours	24 hours	20 hours
Languages	15	15	17	19
Interruption monitor screen	-	-	-	✓
Night mode	-	-	-	✓

For more information or to order, contact your sales representative or visit cardinalhealth.com/KangarooOMNI

© 2023 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, KANGAROO, the Kangaroo LOGO and KANGAROO OMNI are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Patent cardinalhealth.com/patents. Lit. No. 2MS22-1873708 (06/2023)

