



BD[®] Sharps Containers Product catalog

Disposal solutions with a sharp focus
on safety, quality and convenience

● Sharps/Infectious ● Chemotherapy ● Pharmaceutical ● RCRA ● PGII ● BD Recykleen™ ● Accessories





Introducing...

The BD® PGII Sharps Container Portfolio



All-in-one sharps disposal solution

Strict U.S. Department of Transportation (DOT) regulations mandate the use of compliant sharps containers, proper packaging and clear labeling for the safe transport and disposal of Category A and B infectious substances generated in healthcare settings. Packaging for these substances must^{1,2}:

- **Comply** with test standards
- **Employ** a triple-packaging system
- **Display** required UN markings and applicable labeling

What is PGII?

The DOT uses a system known as packing groups (PGs) to classify the level of danger posed by hazardous materials during transport. There are 3 PG classifications³:

- 1 PGI (great danger)
- 2 PGII (medium danger)
- 3 PGIII (minor danger)

BD® PGII Sharps Containers

No boxes. No bags. No bother.



Collect, seal and ship in a single container



Eliminate the additional expense of secondary packaging



Simplify the overall sharps disposal process



Chemotherapy

Pharmaceutical

Hazardous RCRA

Available in 8- and 18-gallon sizes — Sharps/Infectious

1. U.S. Department of Transportation (DOT). Pipeline and Hazardous Materials Safety Administration (PHMSA). Transporting infectious substances safely. Accessed August 8, 2024. <https://www.phmsa.dot.gov/sites/phmsa.dot.gov/files/2024-04/Transporting-Infectious-Substances-Safely-PHH50-0186-0622.pdf> 2. U.S. Department of Transportation (DOT). Pipeline and Hazardous Materials Safety Administration (PHMSA). Package selection marking and labeling requirements. Accessed August 8, 2024. https://www.phmsa.dot.gov/sites/phmsa.dot.gov/files/docs/Package_Selection_Marking_Labeling.pdf 3. U.S. Department of Transportation (DOT) Federal Aviation Administration (FAA). Packaging your dangerous goods. Accessed August 8, 2024. https://www.faa.gov/hazmat/safecargo/how_to_ship/package_for_shipping



Driven by sustainability

We're working smarter to evolve BD Sharps Disposal Solutions

Bringing you high quality products with a sharp focus on a sustainable future

BD is always searching for ways to improve the lifecycle impacts of current and future products. BD Recykleen™ Collectors and Accessories are made from recycled materials.* Our goal is to maximize the recycled content in each of our products. **We are committed to moving the product line in this direction, one step at a time.**

BD is committed to *advancing the world of health™*.



Learn more at
bd.com/sustainability



Evolution in the BD Recykleen™ Collectors brand direction

Products made to suit your clinical needs and the environment

Our goals—high quality and reduced environmental impact

We strive to use as much recycled content as possible in order to deliver high quality disposal solutions with a sharp focus on sustainability.

The BD Recykleen™ Collectors product line offers:

- ✓ Products that are best suited to your clinical needs
- ✓ Products that incorporate recycled content

Counterbalanced door facilitates hands-free disposal

Large opening for easy, safe disposal

VAT locks for secure assembly (visual, audible, tactile)

Handles designed for safe transport

Unique device identification (UDI) compliant



✓ Practical considerations for BD Recykleen™ Collector Products

- **Supply**—availability of recycled resin sufficient to meet the market demand
- **Performance**—same stringent performance requirements for recycled resin as those used for virgin plastic

Contents

○ Innovation and sustainability overview	2
○ Waste stream segregation overview	6
○ BD® Sharps Container Portfolio overview	7
● PGI sharps collectors	8
● Patient room sharps collectors	9
● Large and extra large sharps collectors	10
● Multiuse nestable sharps collectors	11
● Multiuse one-piece sharps collectors	12
● Phlebotomy sharps collectors	13
● Chemotherapy sharps collectors	14
○ Pharmaceutical sharps collectors	15
● RCRA hazardous waste collectors	16
● Trolleys	17
● Cabinets	18
○ Accessories	19
○ BD environmental stewardship	20

Icon key

Look for these symbols throughout the catalog to easily identify key features and benefits:






Waste stream segregation overview*

BD Disposal Solutions offer collectors that meet the disposal needs of the four categories of waste generated in medical settings. There are two main categories of waste: regulated medical waste and hazardous waste. Regulated medical waste includes sharps/infectious, pharmaceutical (non-hazardous) and chemotherapy waste.

The hazardous waste category is for the collection of U.S. Environmental Protection Agency (EPA) Resource Conservation and Recovery Act (RCRA) waste.

This table indicates general guidelines for proper waste stream segregation.



* The graphic offers a standard overview of regulated medical waste disposal. Please refer to applicable local, state and federal regulations for the disposal of regulated medical waste, as each category requires specific handling, transport and treatment to ensure safety and compliance.

† Package or label will indicate.

BD® Sharps Container Portfolio overview



PGII			
	Capacity	Access opening size (inches)	Sample products that fit in collector
Large and extra large	8 gal duo port access	6.4 diameter opening	All items mentioned in the regulated medical waste color categories below, within the large and extra large standard product sizes
	18 gal duo port access	6.4 diameter opening	
Sharps			
Patient room	5.4 qt with counterbalanced door	2.0 × 8.0	Blood collection tubes and collection assembly, syringes up to 60 mL with needle attached
	2 gal with counterbalanced door	2.0 × 8.0	
	3 gal with counterbalanced door	2.0 × 8.0	
Phlebotomy	1 qt	2.5 × 1.2 open (305635)	1 mL–10 mL syringes, blood collection tubes and assembly, slides, pipettes
	1.5 qt	1.5 × 1.5 open (305487)	
Multiuse nestables	8 qt	2.5 × 1.75 diameter open (305343) only	1 mL–20 mL syringes with needle, blood collection tubes, slides, pipettes
	8 qt	1.4 × 1.0 diameter petals (305344, 305059)	
	14 qt	4.22 diameter petals (305456)	
	14 qt	1.8 × 1.1 diameter petals (305464)	
	14 qt	1.0 diameter petals (305480)	
	6 gal	4.22 diameter open (305457)	
6 gal	1.0 diameter petals (305465)		
Multiuse one-piece	1.4 qt	1.8 diameter petals	Small syringes, needles, pipettes, slides
	3.3 qt	1.8 diameter petals	
	6.9 qt	1.8 diameter petals	
	8.2 qt	1.8 diameter petals	
Large and extra large	5 gal	4.2 diameter open	All syringe sizes with needle attached, slides, urine collection assembly, blood culture bottles, blood collection assembly
	9 gal hinge top	9.0 × 10.0	
	9 gal slide top	5.0 × 7.5	
	17 gal hinge top	10.5 × 12.5	
	19 gal slide top	6.0 × 8.5	
Trace chemotherapy			
Multiuse one-piece	3 gal	4.22 diameter open (305076)	All syringe sizes with needle attached, slides, urine collection assembly, blood culture bottles, blood collection assembly
	5 gal	4.2 diameter open (305493)	
Large and extra large	9 gal hinge top	9.0 × 10.0	All syringe sizes, all IV bag sizes with tubing attached, all medication vials, ampules, TPN bags, IV glass bottles
	9 gal slide top	5.0 × 7.5	
	17 gal hinge top	10.5 × 12.5	
	19 gal slide top	6.0 × 8.5	
Non-hazardous pharmaceutical			
Patient room	5.4 qt with counterbalanced door*	2.0 × 8.0	All syringe sizes with needle attached, small IV bags 50 mL–250 mL, medication vials, ampules
	3 gal with counterbalanced door	2.0 × 8.0	All syringe sizes with needle attached, small IV bags 50 mL–250 mL
Multiuse one-piece	3 gal	4.2 diameter open (305633)	All syringe sizes and bags up to 250 mL
Large and extra large	9 gal hinge top	9.0 × 10.0	All syringe sizes, all IV bag sizes with tubing attached, all medication vials, ampules, TPN bags, IV glass bottles
	9 gal slide top	5.0 × 7.5	
	17 gal hinge top	10.5 × 12.5	
	19 gal slide top	6.0 × 8.5	
Hazardous RCRA			
Multiuse one-piece	3 gal*	4.22 diameter open (305076)	Smaller size of solid RCRA hazardous materials plus liquids if contained in a leak-resistant container
	3 gal*	1.0 diameter petal (305075)	
	5 gal	4.2 diameter open (305067)	
Large and extra large	9 gal hinge top*	9.0 × 10.0	Solid RCRA hazardous materials plus liquids if contained in a leak-resistant container
	9 gal slide top*	5.0 × 7.5	
	17 gal hinge top*	10.5 × 12.5	
	19 gal slide top*	6.0 × 8.5	

PGII sharps collectors



All-in-one sharps disposal solution

BD® PGII Sharps Containers simplify the sharps disposal process by eliminating the additional expense of secondary packaging.* Our PGII-rated containers are available in 2 sizes to efficiently manage the 4 types of medical waste generated in healthcare settings: sharps/infectious, chemotherapy, pharmaceutical and hazardous RCRA.

Key product features include:

- Hinge system keeps the lid connected when open
- Duo port lid and cap provide convenient temporary closure
- Leakproof, hermetic, permanent seal in the lid and cap prevents leakage, promotes safe handling and ensures secure final closure
- 16-point perimeter locking lid system adds security against tampering
- Trolley designed for and compatible with all BD® PGII Sharps Containers



Low-temp incinerator



High-temp incinerator



PGII

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/case	Compatible with
304455	PGII 8 gal duo port access lid, red	17.3 x 12.9 x 13.5	9	304463 – PGII Trolley
304456	PGII 8 gal duo port access lid, yellow	17.3 x 12.9 x 13.5	9	
304457	PGII 8 gal duo port access lid, white	17.3 x 12.9 x 13.5	9	
304458	PGII 8 gal duo port access lid, black	17.3 x 12.9 x 13.5	9	
304459	PGII 18 gal duo port access lid, red	17.3 x 12.9 x 24.6	7	
304460	PGII 18 gal duo port access lid, yellow	17.3 x 12.9 x 24.6	7	
304461	PGII 18 gal duo port access lid, white	17.3 x 12.9 x 24.6	7	
304462	PGII 18 gal duo port access lid, black	17.3 x 12.9 x 24.6	7	
304463	PGII Trolley	14.5 x 15.7 x 40.5	1	All 8 and 18 gal (PGII only)

PG: packaging group; PGII: packaging group II.

* BD® PGII Sharps Containers meet U.S. Department of Transportation (DOT) packaging group II requirements for drop (49 CFR 178.603), stacking (49 CFR 178.606), vibration (49 CFR 178.608) and leakage (49 CFR 173.197) testing. BD® PGII Sharps Containers also meet Y performance standards (PGII and PGIII tests) for solid medical waste under UN3291.

Patient room sharps collectors



BD® Patient Room Collectors are available in BD Recykleen™ Collector brand and BD standard brand, and come in a variety of sizes and colors. All BD® Patient Room Collectors have visible fill-level indicators and handles for safe transportation.

Key product features include:

- Compatible with BD locking cabinets and wall brackets to provide security and prevent unwanted access to container contents (see page 18 for more information about locking cabinets)
- Counterbalanced door facilitates hands-free disposal and is designed to reduce the potential of sharps injuries
- Designed to improve tamper resistance by preventing hand access to the contents
- Available in BD Recykleen™ Collectors brand



305517

305053*



Contains recycled content



Hands-free

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/case	Compatible with
305053	3 gal BD Recykleen™ Sharps Collector counterbalanced door, pearl	13.5 x 6 x 15.9	10	305097, 305099, 305409, 305588
305054	3 gal BD Recykleen™ Sharps Collector counterbalanced door, red	13.5 x 6 x 15.9	10	
305517	5.4 qt counterbalanced door, red	12.1 x 4.8 x 10.9	20	305096, 305098, 305447
305551	5.4 qt counterbalanced door, clear	12.1 x 4.8 x 10.9	20	
305426	5.4 qt counterbalanced door, red	12.1 x 4.8 x 10.9	12	
305427	5.4 qt counterbalanced door, clear	12.1 x 4.8 x 10.9	12	
305422	2 gal counterbalanced door, pearl	13.5 x 6 x 11.7	10	305097, 305099, 305409, 305588
305435	2 gal counterbalanced door, red	13.5 x 6 x 11.7	10	
305436	3 gal counterbalanced door, red	13.5 x 6 x 15.9	10	
305443	5.4 qt counterbalanced side entry, red	11.6 x 4.5 x 11.2	20	305447
305444	5.4 qt counterbalanced side entry, pearl	11.6 x 4.5 x 11.2	20	
305425	5.4 qt counterbalanced side entry, pearl	11.6 x 4.5 x 11.2	12	

Large and extra large sharps collectors



BD® Large and Extra Large Collectors are fully compatible with BD Recykleen™ Trolleys to provide stability and mobility (see accessories on page 19 for details).

Key product features include:

- Temporary and final closure features
- Visual, audible, tactile perimeter locks to ensure that the lid is securely attached to the base and to provide security against tampering
- Double-wall handles for safe transport
- Clear tops to ensure visibility of collector fill level
- Easy accessibility for disposal of large bulky waste
- Dual-opening hinge top designed so that the small or the large slide can remain open and the other locked, based on the size of access opening desired



305610



305616



Slide top



Dual open

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/ case	Compatible with
305615	9 gal hinge top	17.7 x 12.0 x 18.6	8	305090
305601	9 gal hinge top with gasket	17.7 x 12.0 x 18.6	8	
305616	9 gal slide top	17.7 x 12.0 x 19	8	305090, 305091
305602	9 gal slide top with gasket	17.7 x 12.0 x 19	8	
305610	17 gal hinge top	20.4 x 14.7 x 23.5	5	305092
305665	17 gal hinge top with gasket	20.4 x 14.7 x 23.5	5	
305609	19 gal slide top	20.4 x 14.7 x 25.9	5	305092, 305093
305666	19 gal slide top with gasket	20.4 x 14.7 x 25.9	5	

Multiuse nestable sharps collectors



BD® Multiuse Nestable Collectors are available in three sizes and can be used with accessories for added security and stability (see accessories on page 19).

Key product features include:

- Clear tops with multiple opening styles to accommodate different sizes of waste and environments
- Handles designed for safe transport
- Funnel-style lids limit accessibility to the collector contents and minimize spillage if the collector tips over
- Open-top lids feature larger diameter openings to accommodate disposal needs in non-patient/public areas
- Available in BD Recykleen™ brand, made from a minimum of 15% recycled materials



305457

8 qt



Regular funnel clear top hinge cap with petals

8 qt



Open clear top hinge cap

14 qt and 6 gal



Large funnel clear top hinge cap with petals

14 qt



Open clear top with tethered cap



305456 with 305485 bracket



Hands-free



Contains recycled content

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/case	Compatible with
305343	8 qt hinge cap	11.3 x 6.8 x 9.8	24	305485, 305018
305344	8 qt hinge cap with petals	11.3 x 6.8 x 9.8	24	
305456	14 qt hinge cap	12.8 x 8.8 x 11.1	20	305485
305464	14 qt hinge cap with needle ports and petals	12.8 x 8.8 x 11.1	20	
305480	14 qt hinge cap with petals	12.8 x 8.8 x 11.1	20	
305457	6 gal hinge cap	12.8 x 8.8 x 17.1	12	
305465	6 gal hinge cap with petals	12.8 x 8.8 x 17.1	12	
305160	6 gal BD Recykleen™, hinge cap	12.8 x 8.8 x 17.1	12	
305059	8 qt BD Recykleen™, hinge cap with needle ports and petals	11.3 x 6.8 x 9.8	24	

Multiuse one-piece sharps collectors



BD® Multiuse One-piece Collectors are available in five sizes and can be used with accessories for added security.

Key product features include:

- Handles designed for safe transport
- Hinge caps designed to provide temporary and final closure features



305557



305488

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/ case	Compatible with
305557	1.4 qt hinge cap with needle ports and petals, red	3.6 x 3.6 x 7.6	36	
305488	3.3 qt hinge cap with petals	9.1 x 5.3 x 5.6	24	305497
305489	6.9 qt hinge cap with petals	9.4 x 5.3 x 11.6	12	
305491	5 gal hinge cap	10.6 x 7.6 x 17.8	8	
305477	5 gal hinge cap with petals	10.6 x 7.6 x 17.8	8	305497
305100	5 gal tethered cap	10.6 x 7.6 x 18.1	8	

Phlebotomy sharps collectors



BD® Phlebotomy Collectors are available in two sizes and can be used with brackets for added stability.

Key product features include:

- Built-in standard safety and pen needle removal ports designed to prevent contact with sharps during disposal
- Clear tops featuring a clear view of the fill level
- Temporary and final closure features for safety and security



305635



305487 with 305970 bracket

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/case	Compatible with
305635	1 qt with needle ports	4.9 x 4.3 x 6.5	60	305970
305487	1.5 qt with needle ports	5.0 x 5.7 x 9.0	36	305970

Chemotherapy sharps collectors



BD® Chemotherapy Collectors feature a variety of lid designs and sizes for disposal of trace chemotherapy sharps waste (see waste stream segregation on page 6). Trace waste is defined by the Centers for Disease Control and Prevention/ National Institute for Occupational Safety and Health (CDC/NIOSH) as “those that contain less than 3% by weight of the original quantity of hazardous drugs.”*

Key product features include:

- Gaskets or liquid-absorbing materials to prevent leakage
- Visual, audible, tactile perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering
- Temporary closure feature and handles designed for safe transport
- Large and extra large collectors that are fully compatible with BD Recykleen™ Trolleys (see page 17)
- Clear lids with multiple opening styles to accommodate different types of waste and environments



305076



305493



305603



305613



Hands-free



Slide top



Dual open



Absorbent



Low-temp incinerator

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/ case	Compatible with
305076	3 gal hinge cap one piece	10.4 x 7.6 x 12.1	12	305497
305493	5 gal hinge cap one piece	10.6 x 7.6 x 17.8	8	305090
305603	9 gal hinge top with gasket	17.7 x 12.0 x 18.6	8	305090, 305091
305604	9 gal slide top with gasket	17.7 x 12.0 x 19	8	305092
305614	17 gal hinge top with gasket	20.4 x 14.7 x 23.5	5	305092, 305093
305613	19 gal slide top with gasket	20.4 x 14.7 x 25.9	5	

* CDC/NIOSH Alert Publication No. 2004-165.

See accessories on page 19 and trolleys on page 17.

Pharmaceutical sharps collectors



BD® Pharmaceutical Collectors feature a variety of lid designs and sizes to collect and dispose of pharmaceutical sharps waste (see waste stream segregation on page 6).

Key product features include:

- Gaskets or liquid-absorbing materials to prevent leakage
- Visual, audible, tactile perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering
- Handles designed for safe transport
- Dual-opening hinge top designed so that the small or the large side can remain open and the other locked, based on the size of access opening desired
- Large and extra large collectors are fully compatible with BD Recykleen™ Trolleys (see page 17)
- 5.4 quart and 3 gallon pharmaceutical collectors with counterbalanced door:
 - Facilitate hands-free disposal and are designed to reduce the potential for sharps injuries
 - Designed to improve tamper resistance by preventing hand access to the contents
 - Compatible with BD® Locking Cabinets to provide security and prevent unwanted access to collector contents
 - Available in BD Recykleen™ brand



305058



305633



305329



305028



Contains recycled content



Hands-free



Slide top



Dual open



Absorbent



Low-temp incinerator

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/ case	Compatible with
305058	5.4 qt BD Recykleen® Sharps Collector counterbalanced door	12.1 x 4.8 x 10.9	20	305096, 305098, 305447
305622	3 gal counterbalanced door	13.5 x 6 x 15.9	10	305097, 305099, 305409, 305588
305633	3 gal hinge cap	10.4 x 7.6 x 12.1	12	305497
305634	9 gal hinge top with gasket	17.7 x 12.0 x 18.6	8	305090
305329	9 gal slide top with gasket	17.7 x 12.0 x 19	8	305090, 305091
305028	17 gal hinge top with gasket	20.4 x 14.7 x 23.5	5	305092
305328	19 gal slide top with gasket	20.4 x 14.7 x 25.9	5	305092, 305093

RCRA hazardous waste collectors



Non-DOT collectors feature a variety of lid designs and sizes and come with a liquid-absorbent material inside the collector for trace waste.

Key product features include:

- BD Recykleen™ RCRA Collectors made from a minimum of 70% recycled materials
- Clear tops with multiple styles to accommodate different types of waste and environments
- Dual-opening hinge top designed so that one side can remain open and the other locked, based on department needs
- Dual-lock cap (3 and 5 gallon) designed to provide temporary closure for safe intermediate transport and secure final closure with two locks
- Visual, audible, tactile perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering
- Handles designed for safe transport
- Temporary closure feature
- Large and extra large collectors that are fully compatible with BD Recykleen™ Trolleys (see page 17)



305069



305071



305066



305075



Contains recycled content



Hands-free



Slide top



Dual open



Absorbent



High-temp incinerator

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/ case	Compatible with
305075	3 gal BD Recykleen™ Sharps Collector hinge cap with petals	10.4 x 7.6 x 12.1	12	
305068	9 gal BD Recykleen™ Sharps Collector hinge top with gasket	17.7 x 12.0 x 18.6	8	305497
305069	9 gal BD Recykleen™ Sharps Collector slide top with gasket	17.7 x 12.0 x 19	8	
305071	17 gal BD Recykleen™ Sharps Collector hinge top with gasket	20.4 x 14.7 x 23.5	5	305909
305072	19 gal BD Recykleen™ Sharps Collector slide top with gasket	20.4 x 14.7 x 25.9	5	305090, 305091
305066	3 gal hinge cap	10.4 x 7.6 x 12.1	12	305497
305067	5 gal hinge cap	10.6 x 7.6 x 17.8	8	305092, 305093

Our full line of trolleys are available in both standard and BD Recykleen™ for the BD® Large and Extra Large Containers. Trolleys provide stability and mobility, available in basic and foot-operated versions for hands-free operation (non-PGII containers only).

Key product features include:

- Designed to minimize the risk of BD® Collectors tipping over
- Wheels that enable the collector to be conveniently placed where it will be most used
- Locking wheels, ensuring the product is mobile only when you want it to be
- The ability to open and close the collector for hands-free disposal (foot-operated trolley with slide top container lids only)

Collectors not included.



305091



305090



305092



305093



304463



Contains recycled content



Hands-free



PGII

BD Recykleen™ Trolleys

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/case	Compatible with
305090	9 gal BD Recykleen™ Basic Trolley	24 x 15.5 x 22	1	All 9 gal
305091	9 gal BD Recykleen™ Foot-Operated Trolley	24 x 15.5 x 25	1	All 9 gal slide top
305092	17 and 19 gal BD Recykleen™ Basic Trolley	26 x 17.4 x 26.6	1	All 17 and 19 gal
305093	19 gal BD Recykleen™ Foot-Operated Trolley	26 x 17.4 x 30	1	All 19 gal slide top

BD® PGII Trolley

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/case	Compatible with
304463	BD® PGII Trolley	14.5 x 15.7 x 40.5	1	All PGII 8 and 18 gal



BD Recykleen™ Cabinets are available in two sizes to secure patient room style containers: small for use with 5.4 quart collectors and large for use with 2 and 3 gallon collectors. Each size is available in both green and brown semi-transparent colors.

Key product features include:

- A universal lock providing security and preventing unwanted access to collectors
- A one-piece transparent door facilitating easy cleaning
- The ability to clearly see the collector fill-level indication

Collectors not included.



Contains recycled content



Security



Hands-free


BD Recykleen™ Cabinets

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/case	Compatible with
305096	5.4 qt BD Recykleen™ Sharps Collector Wall Cabinet, brown	13.3 x 5.6 x 13	1	All 5.4 qt (counterbalanced door only)
305098	5.4 qt BD Recykleen™ Sharps Collector Wall Cabinet, green	13.3 x 5.6 x 13	1	
305097	2 and 3 gal BD Recykleen™ Sharps Collector Wall Cabinet, brown	15 x 7.6 x 18.5	1	All 2/3 gal
305099	2 and 3 gal BD Recykleen™ Sharps Collector Wall Cabinet, green	15 x 7.6 x 18.5	1	



Brackets

Brackets are designed to secure the collector to a wall or cart. Some brackets include additional security via a locking mechanism.



Bracket (305447)



305409



305485

Cabinets

In addition to the BD Recykleen™ brand of cabinets, we offer standard cabinets in a variety of sizes to secure a container to a wall. Many cabinets feature an integrated lock to prevent access to the container contents.



Cabinet (305017)



305970



305018

Brackets and cabinets

Reference number	Product description	Product dimensions (L x W x H inches)	Quantity/case	Compatible with
305475	Locking wall cabinet	7.3 x 6.5 x 11.9	12	3.2 qt
305447	Locking bracket	10.2 x 3.6 x 10.5	10	All 5.4 qt
305409	Locking bracket	13.9 x 4.8 x 4.1	10	2 and 3 gal counterbalanced door
305588	Wire bracket	12.4 x 5.7 x 6	6	2 and 3 gal counterbalanced door
305448	Glove box cabinet bracket	11.9 x 4.5 x 6.3	6	Attaches to cabinets and 5.4 qt locking bracket
305497	Quick release strap	8 x 0.5 x 1.4	12	Multiuse one piece collectors
305017	Locking wall cabinet with shield	13.3 x 5.6 x 13	1	All 5.4 qt (counterbalanced door only)
305550	Wallmate, beige	11.6 x 4.5 x 11.2	1	All 5.4 qt
305420	Locking wall cabinet	15 x 7.6 x 18.5	1	2 and 3 gal counterbalanced door
305451	Large cup bracket	3.6 x 3.7 x 3.5	12	1.4 qt
305018	Wire bracket	10.2 x 7.3 x 8.7	6	8 qt (multi-use nestable)
305494	Padlock	1.0 x 0.4 x 1.7	12	Lockable wall brackets
305970	Small cup bracket	3.6 x 2.9 x 2.8	12	1.0/1.5 qt
305485	Lockable wall bracket	3.5 x 1.4 x 6.0	12	8/14 qt and 6 gal (multi-use nestable)



BD environmental stewardship

- ⊕ The BD approach to sustainability is embedded in our culture and is intrinsic to our purpose, *advancing the world of health™*.

Sustained performance and a dedication to improvement are at our core. Recognizing that human health and a healthy environment are inseparable, BD is committed to reducing our environmental impacts and supporting our customers and suppliers in doing the same.

Learn more at [bd.com](https://www.bd.com)

Contact BD Customer Support at **1.844.8.BD.LIFE** (1.844.823.5433) or your BD Sales Representative for additional information.

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

[bd.com](https://www.bd.com)



BD, the BD Logo, Advancing the world of health and Recykleen are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2025 BD. All rights reserved. BD-2709 (12/25)